



QUALITY CONTROL RANGES

TOTAL GHRELIN RIA KIT CATALOG # GHRT-89HK

STANDARD, LYOPHILIZED CATALOG # 8089-K

STANDARD	TOTAL GHRELIN VALUE
LOT NUMBER: 23T89	6,687 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **2 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # 6089-K

Reconstitute each vial with **1 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	GTTL-123	273 - 566	pg/mL
QC 2	GTTL-223	748 - 1,553	pg/mL

6089-K: GTTL-123.223 (23 September 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA