

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Anti-FAT10 (human specific) antibody produced in rabbit, affinity isolated antibody

Product Number SAB4200836

Product Description

Anti-FAT10 (human specific) antibody is developed in rabbit using a synthetic peptide corresponding to the C-terminal region of human FAT10 (GeneID: 10537), conjugated to KLH as immunogen. The antibody is affinity-purified using the immunizing peptide immobilized on agarose.

Anti-FAT10 antibody specifically recognizes FAT10 from human origin. The antibody may be used in various immunochemical techniques including immunoblotting and immunoprecipitation. Detection of the FAT10 band by Immunoblotting is specifically inhibited by the immunogen.

The Ubiquitin-like protein FAT10, also known as Ubiquitin D (UBD) or Human Leukocyte Antigen (HLA) - F adjacent transcript 10 (FAT10), is a member of the Ubiquitin-like proteins (UBLs) modifier proteins. There are more than 12 identified UBLs including FAT10, SUMOs, NEDD8, ISG15, UFM1, URM1, ATG12, ATG8, FUB1, and HUB1.² UBLs share a similar structure and function, and all UBLs can be covalently attached to their target proteins via the E1–E2–E3 enzymatic cascade in a similar manner to Ubiquitin. UBLylation is a major post-translational modification that in many cases leads to protein degradation in the 26S proteasome.²

FAT10 is a unique UBL containing a C-terminal GG motif available for activation and conjugation. This motif differentiates FAT10 from ubiquitin and other UBLs that need to be cleaved from precursors. Furthermore, unlike ubiquitin that requires a poly-ubiquitin chain, the FAT10 can bind the 26S proteasome as a monomer to mediate proteasomal degradation.²⁻³

Under normal conditions, FAT10 has a short half-life, its expression is repressed by $p53^4$ and is restricted to tissues of the immune system like thymus, lymph nodes, and the spleen. Inflammatory cytokines such as IFN γ and TNF α and downstream STAT3 and NF- κB signaling can induce its expression also in other tissues. $^{2\text{-}3,5}$ FAT10 expression is also upregulated and can be increased up to 100-fold in many cancer types including liver and colon cancer HCC, colorectal, ovarian and uterine carcinomas, triple-negative breast cancer, bladder cancer, gastric cancer, glioma, pancreatic ductal, adenocarcinoma, and osteosarcoma. $^{2,5\text{-}6}$

FATylation has an important role in many cellular process varied from cell maintenance to DNA repair. Understanding FAT10 mechanism and finding its downstream substrates is important to reveal its role in tumorigenesis.¹⁻⁷

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide as a preservative.

Antibody Concentration: ~1.0 mg/mL

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2–8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

<u>Immunoblotting</u>: a working dilution of 1:1,000-1:2,000 is recommended using a whole cell extract of human HepG2 cells treated with TNF α , IFN γ , and MG-132.

Immunoprecipitation: a working concentration of 5-10 μ l/test is recommended using a whole cell extract of human HepG2 cells treated with TNF α , IFN γ , and MG-132.

<u>Note</u>: In order to obtain best results in various techniques and preparations, it is recommended to determine optimal working dilutions by titration test.

References

- 1. Aichem, A., and Groettrup, M., *Int. J. Biochem. Cell Biol.*, **79**, 451-461 (2016).
- 2. Wang, Z., et al., *Mutat. Res.*, **803-805**, 56-75 (2017).
- 3. Schmidtke, G., et al., *Biochim. Biophys. Acta*, **1843**, 97-102 (2014).
- 4. Zhang, D.W., et al., Oncogene, 25, 2318-27 (2006).
- 5. Lukasiak, S., et al., Oncogene, 27, 6068-74 (2008).
- 6. Aichem, A., and Groettrup, M., *Int. J. Biochem. Cell Biol.*, **79**, 451-461 (2016).
- 7. Leng, L., et al., *J. Proteome Res.*, **13**, 268-76 (2014).

VS,DR,OKF,MAM 07/19-1