


Product Information

19511 Mant-ADP triethylammonium salt solution

Molecular Formula:	C ₁₈ H ₂₂ N ₆ O ₁₁ P ₂ (free acid)
Molecular Weight:	560.35 g/mol (free acid)
Form:	clear aqueous solution
Concentration:	10 mM
pH:	7.5 +/-0.5
Spectroscopic Properties:	λ_{\max} 255/355 nm; λ_{exc} 355 nm; λ_{em} 448 nm; ϵ 23.3/5.8 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5)
Package quantity:	150 Units
Unit Definition:	1 unit = 1 ul of a 10 mM solution



Structural formula of Mant-ADP

Applications:MYOSIN-ATPase cycle^[1]SecA cycle^[2]Dissociation kinetic protein kinase A^[3]FRET: kinesin^[4], myosin^[5, 1]Conformational dynamic: DnaC-protein^[6]**Specific Ligands:**Kinesin head domains^[4]Myosin^[5]**Selected References:**

[1] Robertson *et al.* (2005) Structural rearrangements in the active site of smooth muscle myosin. *Biophysical J.* **89**:1882.

[2] Fak *et al.* (2004) Nucleotide exchange from the high-affinity ATP-binding site in SecA is the rate-limiting step in the ATPase cycle of the soluble enzyme and occurs through a specialized conformational state. *Biochemistry* **43**:7307.

[3] Ni *et al.* (2000) Insights into nucleotide binding in protein kinase A using fluorescent adenosine derivatives. *Protein Science* **9**:1818.

[4] Hackney *et al.* (2009) Half-site inhibition of dimeric kinesin head domains by monomeric tail domains. *Biochemistry* **48**:3448.

[5] Sun *et al.* (2006) Dynamics of the upper 50-kDa domain of myosin V examined with fluorescence resonance energy transfer. *J. Biol. Chem.* **281**:5711.

[6] Galletto *et al.* (2005) The nucleotide-binding site of the *Escherichia coli* DnaC protein: Molecular topography of DnaC protein-nucleotide cofactor complex. *Cell Biochem. and Biophys.* **43**:331.

Pinto *et al.* (2011) Structure-activity relationships for the interactions of 2'- and 3'-(O)-(N-methyl)anthraniloyl-substituted purine and pyrimidine nucleotides with mammalian adenylyl cyclases. *Molecular Pharmacology* **82** (4):358.

Chen *et al.* (2009) ADP but Not Pi Dissociation Contributes to Rate Limitation for *Escherichia coli* Rho*. *J. Biol. Chem.* **284** (49):33773.

Taha *et al.* (2009) Molecular Analysis of the Interaction of Anthrax Adenylyl Cyclase Toxin, Edema Factor, with 2'(3')-O-(N-(methyl)anthraniloyl)-Substituted Purine and Pyrimidine Nucleotides. *Molecular Pharmacology* **75** (3):693.

Del Toro Duany *et al.* (2008) The reverse gyrase helicase-like domain is a nucleotide-dependent switch that is attenuated by the topoisomerase domain. *Nucleic Acids Research* **36**(18):5882.

Kainov *et al.* (2008) Structural Basis of Mechanochemical Coupling in a Hexameric Molecular Motor. *The Journal of biological chemistry* **283**(6):3607.

Bujalowski *et al.* (2000) Kinetic mechanism of nucleotide cofactor binding to *Escherichia coli* replicative helicase DnaB protein. stopped-flow kinetic studies using fluorescent, ribose-, and base-modified nucleotide analogues. *Biochemistry* **39**:2106.

Rice *et al.* (1999) A structural change in the kinesin motor protein that drives motility. *Nature* **402**:778.

Cheng *et al.* (1998) Interaction of mant-adenosine nucleotides and magnesium with kinesin. *Biochemistry* **37**:5288.

Bauer *et al.* (1997) X-ray crystal structure and solution fluorescence characterization of Mg-2'(3')-O-(N-methylanthraniloyl) nucleotides bound to the *Dictyostelium discoideum* myosin motor domain. *J. Mol. Biol.* **274**:394.

Bujalowski *et al.* (1994) Structural characteristics of the nucleotide-binding site of *Escherichia coli* primary replicative helicase DnaB protein. Studies with ribose and base-modified fluorescent nucleotide analogs. *Biochemistry* **33**:4682.

Moore *et al.* (1994) Kinetic mechanism of adenine nucleotide binding to and hydrolysis by the *Escherichia coli* Rep monomer. 1. Use of fluorescent nucleotide analogues. *Biochemistry* **33**:14550.

Hiratsuka (1983) New ribose-modified fluorescent analogs of adenine and guanine nucleotides available as substrates for various enzymes. *Biochim. Biophys. Acta* **742**:496.

Storage Conditions: store at -20 °C, protect from light

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

©2016 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.