



**QUALITY CONTROL RANGES
HUMAN GIP (TOTAL) ELISA
CATALOG NO. EZHGIP-54K**

**STANDARD, LYOPHILIZED
CATALOG#E8054-K**

STANDARD	HUMAN GIP VALUE
LOT NUMBER: E54-21K	2076 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

**QUALITY CONTROLS 1 AND 2, LYOPHILIZED
CATALOG #E6054-K**

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HGIP-121	74 - 153	pg/mL
QC 2	HGIP-221	360 - 749	pg/mL

E6054-K: HGIP-121.221 (16 September 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.