



Product Information

SIGMA QUALITY CONTROL TEST PROCEDURE

Enzymatic Assay of UROKINASE (EC 3.4.21.73)

PRINCIPLE:

Plasminogen + H₂O $\xrightarrow{\text{Urokinase}}$ Plasmin

Casein $\xrightarrow{\text{Plasmin}}$ Perchloric Acid Soluble Amino Acids

CONDITIONS: T = 37°C, pH = 7.5, A_{275nm}, Light path = 1 cm

METHOD: Spectrophotometric Stop Rate Determination

REAGENTS:

- A. 60 mM Tris HCl Buffer with 90 mM Sodium Chloride, pH 7.5 at 37°C.
(Prepare 100 ml in deionized water using Trizma Base, Sigma Prod. No. T-1503, and Sodium Chloride, Sigma Prod. No. S-9625. Adjust to pH 7.5 at 37°C with 1 M HCl.)
- B. Porcine Plasminogen Solution (Plas)
(Immediately before use, prepare a solution containing 2 units/ml Plasminogen, Sigma Prod. No. P-1048, in deionized water.)
- C. 1.4% (w/v) α-Casein Suspension (Casein)
(Prepare 100 ml in Reagent A using α-Casein, Sigma Prod. No. C-7891. Do not heat. Stir to make a homogenous suspension.)
- D. 500 mM Perchloric Acid Reagent (PCA)
(Prepare 100 ml in deionized water using Perchloric Acid, Sigma Stock No. 24425-2.)
- E. Urokinase Enzyme Solution
(Immediately before use, prepare a solution containing 1 - 2 units/ml of Urokinase in cold Reagent A.)

**Enzymatic Assay of UROKINASE
(EC 3.4.21.73)**

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable containers:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	0.90	0.90
Reagent B (Plas)	1.00	1.00

Mix by swirling and equilibrate to 37°C. Then add:

Reagent E (Enzyme Solution)	0.10	-----
-----------------------------	------	-------

Immediately following enzyme addition, add:

Reagent C (Casein)	2.00	2.00
--------------------	------	------

Mix by swirling and incubate at 37°C for exactly 15 minutes. Then add:

Reagent D (PCA)	6.00	6.00
Reagent E (Enzyme Solution)	-----	0.10

Mix by swirling and incubate at 25°C for 60 minutes. Filter the solution through a Whatman #50 filter paper¹ and transfer the filtrate to suitable quartz cuvettes. Record the A_{275nm} for both the Test and Blank using a suitable spectrophotometer.

CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(A_{275nm} \text{ Test} - A_{275nm} \text{ Blank})(10)}{(1)(15)(0.1)}$$

10 = Volume (in milliliters) of assay

1 = Change in Absorbance as per the Unit Definition

15 = Time of reaction (in minutes) as per the Unit Definition

0.1 = Volume (in milliliters) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

**Enzymatic Assay of UROKINASE
(EC 3.4.21.73)**

UNIT DEFINITION:

One unit will activate that amount of porcine plasminogen which will produce a $\Delta A_{275\text{nm}}$ of 1.0 per ml per minute at pH 7.5 at 37°C, when measuring perchloric acid soluble products from α -casein (1 cm light path).²

FINAL ASSAY CONCENTRATION:

In a 4 ml reaction mix, the final concentrations are 45 mM Tris, 68 mM sodium chloride, 2 units plasminogen, 0.70% (w/v) α -casein and 0.1 - 0.2 unit urokinase.

REFERENCES:

Hedner, U., Nilsson, I.M., and Robertson, B. (1966) *Thromb. Diath. Haemorrhag.* **16**, 38-50

Lauritsen, O.S. (1966) *Scandinavian Journal of Clinical and Laboratory Investigation* **18**, 73-79

Lauritsen, O.S. (1966) *Scandinavian Journal of Clinical and Laboratory Investigation* **18**, 69-72

Lauritsen, O.S. (1969) *Scandinavian Journal of Clinical and Laboratory Investigation* **23**, 121-128

Lauritsen, O.S. (1968) *Scandinavian Journal of Clinical and Laboratory Investigation* **22**, 239-246

NOTES:

1. Filtering the solution by vacuum or with 0.45 μm syringe filters will result in lower activity.
2. Higher activities are obtained when human plasminogen is used as a substrate.
3. This assay is based on the cited references.
4. All product and stock numbers, unless otherwise indicated, are Sigma product and stock numbers.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.