



QUALITY CONTROL RANGES

MILLIPLEX® MAP Human Neurodegenerative Disease Magnetic Bead Panel 1 Catalog # HNDG1MAG-36K

Control Catalog # HNDG1-6036-1 Control Lot # HNDG1-108.208

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
APOA-1	Control I	HNDG1-108	13 - 28	ng/mL
	Control II	HNDG1-208	66 - 137	ng/mL
APOC-3	Control I	HNDG1-108	0.22 - 0.46	ng/mL
	Control II	HNDG1-208	0.98 - 2.0	ng/mL
APOE	Control I	HNDG1-108	0.93 - 1.7	ng/mL
	Control II	HNDG1-208	4.4 - 9.1	ng/mL
HCFH	Control I	HNDG1-108	6.2 - 13	ng/mL
	Control II	HNDG1-208	31 - 63	ng/mL
HCC3	Control I	HNDG1-108	0.91 - 1.9	ng/mL
	Control II	HNDG1-208	4.3 - 9.0	ng/mL
HA2MG	Control I	HNDG1-108	10.1 - 21	ng/mL
	Control II	HNDG1-208	46 - 96	ng/mL

HNDG1-6036-1: HNDG1-108.208 (20 October 2022)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA