

QUALITY CONTROL RANGES

Milliplex Rat Adipokine Catalog # RADPK-81K

Milliplex Rat Adiponectin Single Plex Catalog # RADPK-81K-ADPN

Milliplex Rat Adipocyte Catalog # RADPCYT-82K

Control Catalog # LRA-6081

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
IL-1 β	Control I	RADIP106	197 – 408	pg/mL
	Control II	RADIP206	763 – 1584	pg/mL
IL-6	Control I	RADIP106	1772 – 3679	pg/mL
	Control II	RADIP206	6227 – 12933	pg/mL
Adiponectin	Control I	RADIP106	9887 – 20535	pg/mL
	Control II	RADIP206	40582 - 84286	pg/mL
Leptin	Control I	RADIP106	377 – 783	pg/mL
	Control II	RADIP206	1295 - 2689	pg/mL
MCP-1	Control I	RADIP106	374 – 777	pg/mL
	Control II	RADIP206	1363 - 2830	pg/mL
Insulin	Control I	RADIP106	781 – 1622	pg/mL
	Control II	RADIP206	2610 - 5421	pg/mL
TNF- α	Control I	RADIP106	52 – 108	pg/mL
	Control II	RADIP206	186 – 401	pg/mL
Total PAI-1	Control I	RADIP106	2580 – 5358	pg/mL
	Control II	RADIP206	8021 – 16660	pg/mL
Active PAI-1	Control I	RADIP106	2371 – 4925	pg/mL
	Control II	RADIP206	9109 - 18919	pg/mL

EMD Millipore • (636) 441-8400 • (636) 441-8050 (FAX) • www.emdmillipore.com