



## QUALITY CONTROL RANGES

### HUMAN AMYLIN ELISA CATALOG # EZHA-52K STANDARD CATALOG #E8051-K

STANDARD Lot Number. 7905-25K	AMYLIN VALUE 163 pM
----------------------------------	------------------------

The exact value of lyophilized standard is shown above.

Reconstitute with **1.0 mL** distilled or deionized water.

### QUALITY CONTROLS 1 AND 2 CATALOG #E6051-K

Reconstitute each vial with **0.25 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	EA-125	10 - 20	pM
QC 2	EA-225	48 - 99	pM

E6051-K: EA-125.225 (22 AUG 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • SigmaAldrich.com

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.