



Saint Louis, Missouri 63103 USA  
Telephone (800) 325-5832 (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## ProductInformation

### ANTI-CAT WHOLE SERUM

Developed in Rabbit

Delipidized, Whole Antiserum

Product No. **C7415**

#### Product Description

Anti-Cat Whole Serum was developed in rabbit using pooled normal cat serum as the immunogen. The anti-serum has been treated to remove lipoproteins.

Strong reactivity with normal cat serum has been determined by immunoelectrophoresis (IEP). This antiserum has not been assayed for interspecies cross reactivity.

#### Reagents

Rabbit anti-Cat Serum is provided as a liquid containing 0.1% sodium azide (see MSDS)\* as preservative.

#### Precautions

\*Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

#### Storage

For continuous use, store at 2-8EC for a maximum of one month. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Pcs12/99

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.