

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

## **ProductInformation**

## ACETAMINOPHEN SIGMA REFERENCE STANDARD

Product No. A3035

Lot No: 46K8802 CAS No: 103-90-2

Storage Temperature: room temperature

Synonyms: Acetamide, *N*-(4-hydroxyphenyl)-.

4'-Hydroxyacetanilide.

Formula: C<sub>8</sub>H<sub>9</sub>NO<sub>2</sub>

Formula Weight: 151.16

#### **Product Description**

This product meets USP specifications (Current through USP 29, Supplement 1) and is traceable to USP Reference Standard lot J2C423.


Store vial at room temperature. Protect from light. Dry a portion over silica gel for 18 hours before using. Use promptly. Discard unused material.


TEST	USP 29 SPECIFICATIONS	Result
		Sigma R.S. Lot 46K8802
Infrared absorption spectrum KBr dispersion (197K)	compares to standard	conforms to USP reference standard lot J2C423
Ultraviolet absorption spectrum 5 μg/mL 0.1 N HCl in methanol ⟨197U⟩	compares to standard	conforms to USP reference standard lot J2C423
Identification C	meets requirement	passes
Melting range ⟨741⟩	between 168° and 172°	168.7°-168.9°
Water, Method I (921)	not more than 0.5%	passes
Residue on ignition (281)	not more than 0.1%	passes
Chloride ⟨221⟩	not more than 0.014%	passes
Sulfate (221)	not more than 0.02%	passes
Sulfide	no coloration or spotting of test paper	passes
Heavy metals, Method II (231)	not more than 0.001%	passes
Free <i>p</i> -aminophenol	not more than 0.005%	passes
Limit of <i>p</i> -chloroacetanilide	not more than 0.001%	passes
Readily carbonizable substances (271)	the test solution has no more color than <i>Matching Fluid A</i>	passes
Assay	not less than 98.0 percent and not more than 101.0 percent C <sub>8</sub> H <sub>9</sub> NO <sub>2</sub> , calculated on the anhydrous basis	100.2%

See reverse side for IR & UV spectra of Sigma Reference Standard.

## ACETAMINOPHEN SIGMA REFERENCE STANDARD


Product No. **A3035** Lot 46K8802





# ULTRAVIOLET SPECTRUM

#### $5.16~\mu\text{g/mL}$ in 0.1 N HCl in methanol



Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.