



## Product Information

### IL-2/Fc CHIMERA (NON-LYTIC)

Mouse, Recombinant  
Expressed in CHO cells

Product Number **I 9905**

Storage Temperature  $-20\text{ }^{\circ}\text{C}$

### Product Description

IL-2 (Interleukin-2)/Fc Chimera is a soluble 94 kDa dimeric fusion protein consisting of mouse IL-2 fused to mutant mouse Fc $\gamma$ 2a Fc. It is purified from tissue culture supernatants of CHO cell transfectants. This fusion protein possesses both the biological functions of IL-2 as an immune anti-inflammatory agent and the prolonged circulating half-life determined by the Fc domain. Mutations to the complement (C1q) and Fc $\gamma$ RI binding sites of the Fc $\gamma$ 2a fragment render IL-10/Fc unable to facilitate antibody directed cytotoxicity (ADCC) and complement mediated cytotoxicity (CDC).<sup>1,2</sup>

Interleukin-2 (IL-2), also known as T cell growth factor, is a 133 amino acid glycoprotein that is secreted by activated T lymphocytes and induces proliferation and maturation of activated T cells, natural killer cells, lymphokine activated killer cells, monocytes, macrophages, and oligodendrocytes.<sup>3</sup> IL-2 stimulates the proliferation of antibody-producing B cells, activates neutrophils, modulates the expression of the IL-2 receptor, and induces mononuclear cells to secrete  $\gamma$ -interferon and TNF- $\alpha$  and  $-\beta$ .<sup>3-6</sup> In addition, recent studies have shown that IL-2 is required for activation-induced apoptosis, an important homeostasis mechanism for the maintenance of peripheral tolerance to self antigens.<sup>7,8</sup>

### Reagent

IL-2/Fc chimera is supplied as 10  $\mu\text{g}$  of protein in 0.22  $\mu\text{m}$  sterile-filtered PBS, pH 7.4 (50 mM sodium phosphate, 100 mM potassium chloride, 150 mM sodium chloride) and contains no preservatives.

### Preparation Instructions

The chimera solution can be further diluted to the desired working concentration in sterile PBS or culture medium.

### Storage/Stability

Store at  $-20\text{ }^{\circ}\text{C}$ . Store working solutions at  $4\text{ }^{\circ}\text{C}$  for up to one week. Repeated freeze/thaw cycles are not recommended.

### Product Profile

The biological activity of recombinant IL-2/Fc is determined in a proliferation assay using CTLL-2 indicator cells.<sup>9</sup>

Specific Activity: 1 to  $1.5 \times 10^6$  Units/mg

A Unit is defined using rIL-2 as the reference in the CTLL-2 cell proliferation assay.

Purity: >98% by SDS-PAGE

### References

1. Zheng, X.X., et al., J. Immunol., **163**, 4041 (1999).
2. Zheng, X.X., et al., J. Immunol., **154**, 5590 (1995).
3. Smith, K., Science, **240**, 1169 (1988).
4. Robb, R.J., et al., Proc. Natl. Acad. Sci. USA **81**, 6486 (1984).
5. Gillis, S., et al., J. Immunol. **120**, 2027 (1978).
6. Fletcher, M., and Goldstein, A.L., Lymphokine Res. **6**, 45 (1987).
7. Lenardo, K.M., Nature **353**, 858 (1991).
8. Kneitz, B., et al., Eur. J. Immunol. **25**, 2572 (1995).
9. Gearing, A. et al., In: Lymphokines and Interferons, A Practical Approach, Clemens, M. et al. (eds.), p. 296 (IRL Press, Oxford, 1987).

CMH 9/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.