



VQUALITY CONTROL RANGES

PIIANP ELISA CATALOG # EZPIIANP-53K

STANDARD, LYOPHILIZED CATALOG#E8053-K

| | |
|---------------------|--------------|
| STANDARD | PIIANP VALUE |
| Lot Number: E53-13K | 3295 ng/mL |

The exact value of lyophilized standard is shown above.

Reconstitute with **0.2 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED

Reconstitute each vial with **0.2 mL** distilled or deionized water.

| CONTROL | LOT NUMBER | EXPECTED RANGE | UNITS |
|---------|------------|----------------|-------|
| QC 1 | PIIA-113 | 85 - 176 | ng/mL |
| QC 2 | PIIA-213 | 393 - 1170 | ng/mL |

E6053-K: PIIA-113.213 (02 October 2020)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.