

QUALITY CONTROL RANGES
Human Adipocyte Magnetic Bead Panel
Catalog # HADCYMAG-61K
QC Catalog #: HADCYT-6061



Control 1 Lot#: HADCY-107A / Control 2 Lot#: HADCY-207A

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
NGF	Control 1	23 - 49	pg/mL
	Control 2	122 - 253	pg/mL
IL-6	Control 1	67 - 140	pg/mL
	Control 2	372 - 773	pg/mL
LEPTIN	Control 1	832 - 1729	pg/mL
	Control 2	4104 - 8525	pg/mL
IL-8	Control 1	19 - 39	pg/mL
	Control 2	99 - 206	pg/mL
HGF	Control 1	502 - 1043	pg/mL
	Control 2	2247 - 4668	pg/mL
ADPN	Control 1	1765 - 3666	pg/mL
	Control 2	9776 - 20304	pg/mL
MCP-1	Control 1	80 - 167	pg/mL
	Control 2	366 - 759	pg/mL
TNF- α	Control 1	28 - 59	pg/mL
	Control 2	147 - 305	pg/mL
RESISTIN	Control 1	463 - 962	pg/mL
	Control 2	2438 - 5064	pg/mL
IL-1B	Control 1	38 - 79	pg/mL
	Control 2	193 - 401	pg/mL
PAI-1(total)	Control 1	734 - 1525	pg/mL
	Control 2	3825 - 7945	pg/mL

HADCYT-6061: HADCY-107A.207A (02 August 2021)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.