

Product Information Sheet

Interleukin-18 from rat

recombinant, expressed in *E. coli***I0531**Storage Temperature $-20\text{ }^{\circ}\text{C}$

Product Description

Interleukin-18 (IL-18) shares some biological activities with IL-12 and structural similarities with the IL-1 family. Structurally, IL-18 and IL-1 β share primary amino acids sequences (signature sequence motif) and are similarly folded as all- β pleated sheet molecules. Similar to IL-1 β , production of IL-18 requires caspase-1 (interleukin-1 β converting enzyme, ICE) to cleave the pro-domain from the precursor protein to produce active mature IL-18.

The primary immune actions of IL-18 appear to be through the induction of high levels of IFN- γ production. In the T cell helper type 1 (Th1) response, IL-18 induces IFN- γ production in T cells and NK cells. Gene expression and synthesis of TNF, IL-1, Fas ligand, and several chemokines are induced by IL-18. IL-18 selectively enhances the FasL-mediated cytotoxicity of Th1, but not Th0 or Th2 cells. IL-18 induces activated B cells to produce IFN- γ that inhibits IgE production. It also plays a key role in immunoregulation and participates in both innate and acquired immunity.

IL-18 has been detected in keratinocytes, small intestine epithelial cells, adrenal cells, macrophages, pancreas, skeletal muscle, liver, lung, and PBMCs (peripheral blood mononuclear cells). IL-18 actions are mediated by an IL-18 receptor (IL-18 R) complex comprised of a binding chain (IL-18 R α , also called IL-1 receptor-related protein, IL-1 Rrp) and a signal-transducing chain.

This recombinant Rat Interleukin-18 (IL-18) product is produced from a DNA sequence encoding the mature rat IL-18 protein. It is supplied lyophilized from a 0.2 μm filtered solution of 20 mM MOPS, 50 mM sodium sulfate, 0.5 mM EDTA, and 0.5 mM DTT, pH 7.2, containing 50 μg bovine serum albumin per 1 μg of cytokine with 5% trehalose.

The methionyl form of recombinant rat IL-18, a 159 amino acid protein, has a predicted molecular mass of ~18.4 kDa.

Purity

> 97% (SDS-PAGE)

ED₅₀: typically 10–60 ng/mL

Biological activity is measured by the ability to induce mouse IFN- γ production by activated mouse T cells. The ED₅₀ is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell-based bioassay.

Endotoxin level

<0.1 EU/ μg of the cytokine

(LAL [Limulus amoebocyte lysate] method)

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The product ships on dry ice and storage at $-20\text{ }^{\circ}\text{C}$ is recommended. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Do not store in a frost-free freezer.

References

1. Conti, B. et al., Induction of interferon- γ -inducing factor in the adrenal cortex. *J. Biol. Chem.*, 272, 2035 (1997).
2. Dinarello, C.A., Interleukin-18. *Methods*, 19, 121 (1999).

3. Antonysamy, M.A. et al., Interleukins 16, 17 and 18, in *The Cytokine Handbook*, 3rd Edition, Thomson, A.W., ed., pp. 465-489 (Academic Press, San Diego, 1998).
4. Bazan, J.F. et al., A newly defined interleukin-1 *Nature*, 379, 591 (1996).
5. Kohno, K. et al., IFN-gamma-inducing factor (IGIF) is a costimulatory factor on the activation of Th1 but not Th2 cells and exerts its effect independently of IL-12. *J. Immunol.*, 158, 1541 (1997).
6. Gu, Y. et al., Activation of interferon- γ inducing factor mediated by interleukin-1 β converting enzyme. *Science*, 275, 206 (1997).
7. Ghayur, T. et al., Caspase-1 processes IFN- γ -inducing factor and regulates LPS-induced IFN- γ production. *Nature*, 386, 619 (1997).
8. Micallef, M.J. et al., Interferon- γ -inducing factor enhances T helper 1 cytokine production by stimulated human T cells: synergism with interleukin-12 for interferon- γ production. *Eur. J. Immunol.*, 26, 1647 (1996).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

I0531pis Rev 03/22

**Millipore
Sigma**