

50195 Listeria UVM I Supplement (Listeria primary selective enrichment supplement UVM I)

An antimicrobial supplement recommended for selective isolation and cultivation of Listeria species.

Composition

(per vial, sufficient for 500 ml medium)

Acridine hydrochloride	6.0 mg
Nalidixic acid	10.0 mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, cooled (45-50 °C) UVM Listeria Selective Enrichment Broth, modified I & II. (Cat. No. 90554). Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

