

## shelf life information data sheet

**Product Family: Polygard® CE**

### shelf life information

Shelf Life: 60 Months

Expiration Date: N/A

Storage Conditions: Ambient temperature

Applicable Catalog Number(s): See Attached list

### quality attributes studied

Visual Examination

Water Flow Rate

Not all products listed have undergone stability testing to determine the products' shelf lives. The design of the stability schedule is such that only select product(s) within a product family are tested. The design assumes that the products selected are representative of the family as they utilize the same materials of construction, production processes, and packaging components as the products within the family.

The product(s) above does not have a published expiration date or shelf life claim. The shelf life stated above is a guideline based upon available product stability data which is constantly updated as new data is made available. It is not a guarantee or a claim designating the time prior to which a lot of the product will remain within the shelf life specification if stored under the defined conditions, and after which it must not be put into use. It is not an expiration date or a use-by-date. Further, as this is not a product claim, updates to product shelf lives do not require customer notification.

For firms that must establish expiration dates or use-by-dates for this product(s) utilized in a cGMP manufacturing environment, it is recommended that the firms assign a default use-by-date based on their internal procedures/practices.

A handwritten signature in black ink, appearing to read "T. Brackett".

28-Oct-19

Thomas Brackett  
Supplier Quality Engineering Manager  
Third Party Finished Goods Quality

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.



---

CEA571S06
CEA572S06
CEA573S06
CEA574S06
CEA501S06
CEA502S06
CEA503S06
CEA504S06
CEA5F1S06
CEA5F2S06
CEA5F3S06
CEA5F4S06
CE0171S06
CE0172S06
CE0173S06
CE0174S06
CE0101S06
CE0102S06
CE0103S06
CE0104S06
CE01F1S06
CE01F2S06
CE01F3S06
CE01F4S06
CE0271S06
CE0272S06
CE0273S06
CE0274S06
CE0201S06
CE0202S06
CE0203S06
CE0204S06
CE02F1S06
CE02F2S06
CE02F3S06
CE02F4S06