



# ORANGE IS THE NEW BLUE

Avoid health risks that occur when working with toxic and carcinogenic blue silica gel.

Based on iron salt, our orange silica gels give you nothing but reliable results. Explore our range of orange silica gels and many other non-toxic drying agents – to enjoy safety and quality perfected for your intended use.



Explore our complete range  
of drying agents on:  
**[SigmaAldrich.com/drying-agents](https://SigmaAldrich.com/drying-agents)**




The life science  
business of Merck  
operates as  
MilliporeSigma in  
the U.S. and Canada.

**Supelco®**  
Analytical Products

# Ordering information

| Product                                                                          | Content                                                             | Packaging           | Ord. No.            |
|----------------------------------------------------------------------------------|---------------------------------------------------------------------|---------------------|---------------------|
| <b>Silica gel [SiO<sub>2</sub>]</b>                                              |                                                                     |                     |                     |
|  | Silica gel granules, desiccant ~ 0.2-1 mm                           | 1 kg Plastic bottle | <b>1.01905.1000</b> |
|                                                                                  | Silica gel granules, desiccant ~ 2-5 mm                             | 1 kg Plastic bottle | <b>1.01907.1000</b> |
|                                                                                  |                                                                     | 5 kg Plastic bottle | <b>1.01907.5000</b> |
|                                                                                  | Silica gel with moisture indicator (brown gel), desiccant ~ 1-4 mm  | 1 kg Plastic bottle | <b>1.01972.1000</b> |
|                                                                                  |                                                                     | 5 kg Plastic bottle | <b>1.01972.5000</b> |
|                                                                                  |                                                                     | 25 kg Plastic drum  | <b>1.01972.9025</b> |
|                                                                                  | Silica gel with moisture indicator (orange gel), granulate ~ 1-3 mm | 1 kg Plastic bottle | <b>1.01969.1000</b> |
|                                                                                  |                                                                     | 5 kg Plastic bottle | <b>1.01969.5000</b> |
|                                                                                  |                                                                     | 25 kg Plastic drum  | <b>1.01969.9025</b> |
|                                                                                  | Silica gel beads, desiccant ~ 2-5 mm                                | 1 kg Plastic bottle | <b>1.07735.1000</b> |

|                                                                                  |                                                                                              |                         |                     |
|----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|-------------------------|---------------------|
| <b>Desiccant sachets [SiO<sub>2</sub>]</b>                                       |                                                                                              |                         |                     |
|  | Desiccant sachet 3 g silica gel with humidity indicator (orange gel), sachet: 4 x 7 cm       | 100 units Metal can     | <b>1.03803.0001</b> |
|                                                                                  |                                                                                              | 1000 units Fiber carton | <b>1.03803.0002</b> |
|                                                                                  | Desiccant sachet 10 g silica gel with humidity indicator (orange gel), sachet: 7 x 9 cm      | 50 units Metal can      | <b>1.03804.0001</b> |
|                                                                                  | Desiccant sachet 100 g silica gel with humidity indicator (orange gel), sachet: 15 x 14 cm   | 10 units Metal can      | <b>1.03805.0001</b> |
|                                                                                  | Desiccant sachet 250 g silica gel with humidity indicator (orange gel), sachet: 15 x 20.5 cm | 10 units Metal can      | <b>1.03806.0001</b> |

|                                                                                                                                                                                                                                                                   |                                                                                                           |                      |                     |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|----------------------|---------------------|
| <b>Molecular sieves</b>                                                                                                                                                                                                                                           |                                                                                                           |                      |                     |
| <br><br> | Molecular sieve 0.3 nm beads ~ 2 mm (suitable for use in Karl Fischer titration)                          | 250 g Plastic bottle | <b>1.05704.0250</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 1 kg Plastic bottle  | <b>1.05704.1000</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 10 kg Plastic bucket | <b>1.05704.9010</b> |
|                                                                                                                                                                                                                                                                   | Molecular sieve 0.4 nm beads ~ 2 mm Reag. Ph Eur                                                          | 250 g Plastic bottle | <b>1.05708.0250</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 1 kg Plastic bottle  | <b>1.05708.1000</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 10 kg Plastic bucket | <b>1.05708.9010</b> |
|                                                                                                                                                                                                                                                                   | Molecular sieve 1.0 nm beads ~ 2 mm                                                                       | 1 kg Glass bottle    | <b>1.05703.1000</b> |
|                                                                                                                                                                                                                                                                   | Molecular sieve 0.3 nm beads, with moisture indicator ~ 2 mm (suitable for use in Karl Fischer titration) | 250 g Plastic bottle | <b>1.05734.0250</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 1 kg Plastic bottle  | <b>1.05734.1000</b> |
|                                                                                                                                                                                                                                                                   | Molecular sieve 0.4 nm beads, with moisture indicator ~ 2 mm                                              | 250 g Glass bottle   | <b>1.05739.0250</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 1 kg Glass bottle    | <b>1.05739.1000</b> |
|                                                                                                                                                                                                                                                                   | Molecular sieve 0.3 nm rods ~ 1.6 mm (1/16")                                                              | 250 g Plastic bottle | <b>1.05741.0250</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 1 kg Plastic bottle  | <b>1.05741.1000</b> |
|                                                                                                                                                                                                                                                                   |                                                                                                           | 10 kg Plastic bucket | <b>1.05741.9010</b> |
|                                                                                                                                                                                                                                                                   | Molecular sieve 0.4 nm rods ~ 1.6 mm (1/16")                                                              | 1 kg Plastic bottle  | <b>1.05743.1000</b> |

The smoothest way to seal & lubricate your desiccator:



our **silicone grease**

Art.No. 1.07746

© 2019 Merck KGaA, Darmstadt, Germany and its affiliates. All Rights Reserved. Merck, the vibrant M, and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MK\_FL4274EN Ver. 1.0  
07/2019

