



## QUALITY CONTROL RANGES

### RAT / MOUSE GHRELIN (TOTAL) ELISA KIT CATALOG # EZRGRT-91K

#### STANDARD, LYOPHILIZED CATALOG # E8091-K

STANDARD	TOTAL GHRELIN VALUE
Lot Number: ERGA-15K	10 ng/mL

Hydrate standard with **2.0mL** deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6091-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	RAGR-115	0.29 – 0.61	ng/mL
QC 2	RAGR-215	2.6 – 5.5	ng/mL

E6091-K: RAGR-115.215 (27 December 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA