



## RABBIT ANTI-DYNAMITIN (p50), aa 1-83 POLYCLONAL ANTIBODY

<b>CATALOG NUMBER:</b>	AB5869
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	50 µL
<b>SPECIFICITY:</b>	Dynamitin (p50), aa 1-83
<b>IMMUNOGEN:</b>	Synthetic peptide from human Dynamitin, aa 1-83
<b>APPLICATIONS:</b>	Western blot: 1:1,000 Optimal working dilutions must be determined by the end user.
<b>SPECIES REACTIVITY:</b>	Human, rat, mouse and rabbit. Other species have not been tested.
<b>FORMAT:</b>	Purified immunoglobulin.
<b>PRESENTATION:</b>	Liquid in PBS with 50% Glycerol and 0.02% sodium azide.
<b>STORAGE/HANDLING:</b>	Maintain at -20°C in undiluted for up to 6 months from date of receipt. Avoid repeated freeze/thaw cycles.

*For research use only; not for use as a diagnostic.*

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*