

XXXXXX

7.76071.0011 03/2025

1.16995

Supelco®

sigmaaldrich.com

50 Tests


Reflectoquant®

**Nitrate Test** for RQflex® 20**Nitrat-Test** für RQflex® 203 - 90 mg/l NO<sub>3</sub><sup>-</sup>

Reaction time: 60 sec



Warning



Please consult Instruction for Use for further information and observe the RQflex® manual.



For IfU please scan code with My M Safety App or go to [sigmaaldrich.com](http://sigmaaldrich.com)



©2025

Store cold and dry • (2 - 8 °C) • Kalt und trocken lagern

Lot

min. shelf life (YYYY/MM/DD)

Merck KGaA, 64271 Darmstadt, Germany  
EMD Millipore Corporation, Burlington MA 01803, USA  
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,  
Ontario, L6H 6J8, Canada