



**QUALITY CONTROL RANGES**  
**Milliplex Canine Kidney Toxicity Magnetic Bead Panel 2**  
**Catalog # CKT2MAG-97K**  
**Control Catalog # CKT2-6097-2**

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Albumin	Control I	CKT2-107	214 – 444	ng/mL
	Control II	CKT2-207	2022 - 4199	ng/mL
b-2-Microglobulin	Control I	CKT2-107	0.34 - 0.70	ng/mL
	Control II	CKT2-207	3.3 - 7.0	ng/mL
TFF-3	Control I	CKT2-107	0.37 - 0.78	ng/mL
	Control II	CKT2-207	2.9 - 6.1	ng/mL
RBP4	Control I	CKT2-107	0.17 - 0.35	ng/mL
	Control II	CKT2-207	1.7 - 3.5	ng/mL

**Note: The Quality Control Ranges are generated using the same day assay format. Quality Control values in culture media are not tested.**

CKT2-6097-2: CKT2-107.207 (12 October 2022)

EMD MILIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA