



SIGMA QUALITY CONTROL TEST

Product Information

Enzymatic Assay of NUCLEOSIDE PHOSPHORYLASE (EC 2.4.2.1)

PRINCIPLE:

Inosine + Pi $\xrightarrow{\text{Nucleoside Phosphorylase}}$ Hypoxanthine + Ribose-1-PO₄

Hypoxanthine + 2H₂O + 2O₂ $\xrightarrow{\text{Xanthine Oxidase}}$ Uric Acid + 2H₂O₂

Abbreviations:

Pi = Inorganic Phosphate

CONDITIONS: T = 25°C, pH = 7.4, A_{293nm}, Light path = 1 cm

METHOD: Continuous Spectrophotometric Rate Determination

REAGENTS:

- A. 100 mM Potassium Phosphate Buffer, pH 7.4 at 25°C
(Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Prod. No. P-5379. Adjust to pH 7.4 at 25°C with 1 M NaOH.)
- B. 7.5 mM Inosine Solution
(Prepare 5 ml in deionized water using Inosine, Prod. No. I-4125.)
- C. Xanthine Oxidase Enzyme Solution
(Immediately before use, prepare a solution containing 10 units/ml of Xanthine Oxidase, Prod. No. X-4500, in ice cold deionized water.)
- D. Nucleoside Phosphorylase Enzyme Solution
(Immediately before use, prepare a solution containing 0.125 units/ml of Nucleoside Phosphorylase in ice cold deionized water.)

**Enzymatic Assay of NUCLEOSIDE PHOSPHORYLASE
(EC 2.4.2.1)**

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	2.70	2.70
Reagent B (Inosine)	0.10	0.10
Reagent C (Xanthine Oxidase)	0.10	0.10

Mix by inversion and equilibrate to 25°C. Monitor the $A_{293\text{nm}}$ until constant, using a suitably thermostatted spectrophotometer. Then add:

Reagent D (Nucleoside Phosphorylase)	0.10	-----
Deionized Water	-----	0.10

Immediately mix by inversion and record the increase in $A_{293\text{nm}}$ for approximately 5 minutes. Obtain the $\Delta A_{293\text{nm}}$ /minute using the maximum linear rate for both the Test and Blank.

CALCULATIONS:

$$\text{Units/mg enzyme} = \frac{(\Delta A_{293\text{nm}}/\text{min Test} - \Delta A_{293\text{nm}}/\text{min Blank})}{(12.0) (\text{mg enzyme/ml RM})}$$

12.0 = Millimolar extinction coefficient of Uric Acid at 293nm
RM = Reaction Mix

UNIT DEFINITION:

One unit will cause the phosphorolysis of 1.0 μmole of inosine to hypoxanthine and ribose 1-phosphate per minute at pH 7.4 at 25°C.

FINAL ASSAY CONCENTRATION:

In a 3.0 ml reaction mix, the final concentrations are 90 mM potassium phosphate, 0.25 mM inosine, 1.0 units of xanthine oxidase and 0.0125 units of nucleoside phosphorylase.

**Enzymatic Assay of NUCLEOSIDE PHOSPHORYLASE
(EC 2.4.2.1)**

NOTES:

1. Xanthine Oxidase - One unit will convert 1.0 μ mole of xanthine to uric acid per minute at pH 7.5 at 25°C.
2. All products and stock numbers, unless otherwise indicated, are Sigma product and stock numbers.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.