



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Anti-SMC1L1 (C-terminal)

produced in rabbit, IgG fraction of antiserum

Catalog Number **S6446**

Product Description

Anti SMC1L1 (C-terminal) is produced in rabbit using a synthetic peptide corresponding to amino acids 1210-1225 of human SMC1L1, conjugated to KLH via an N-terminal cysteine residue, as immunogen. The immunizing sequence is conserved in mouse and rat. Whole antiserum is fractionated and then further purified by ion-exchange chromatography to provide the IgG fraction of antiserum that is essentially free of other rabbit serum proteins.

Anti SMC1L1 (also known as SMC1 and SMC1 α) specifically recognizes SMC1L1 by immunoblotting (160 kDa), immunoprecipitation, and immunofluorescence. Staining of the SMC1L1 bands in immunoblotting are specifically inhibited by the immunizing peptide.

Proper cohesion of sister chromatids is a prerequisite for the correct segregation of chromosomes during cell division. The cohesin chromosome complex is required for sister chromatid cohesion.¹ This complex is partly composed of two Structural Maintenance of Chromosomes (SMC) proteins, SMC3, and either SMC1L1 (SMC1, SMC1 α) or SMC1L2. The diverse functions of the SMC1 proteins range far beyond chromosome segregation and involve nearly all aspects of chromosome behavior including DNA repair.¹⁻³ SMC1 and SMC2 were originally identified in *Saccharomyces cerevisiae* as genes required for proper condensation and segregation of chromosomes.⁴ There are at least six SMC family members that form three heterodimers in specific combinations: SMC1L1 (SMC1) and SMC3 constitute the core of the cohesin complex that maintains sister chromatid cohesion, whereas SMC2 and SMC4 are components of the condensin complexes that mediate chromosome condensation during mitosis. A third complex, containing SMC5 and SMC6, is mainly involved in the cellular response to DNA damage in yeast as well as humans.^{1,5,6} Although SMC1L1 was originally identified as a component of cohesin, it has

been shown to be phosphorylated by ATM indicating a potential role for this protein in DNA repair. ATM (ataxia telangiectasia mutated) and ATR (ataxia telangiectasia and Rad3 related) are key signal transducers of the DNA damage response. SMC1L1 is phosphorylated at Ser⁹⁵⁷ and Ser⁹⁶⁶ in an ATM-dependent manner after ionizing irradiation in response to DNA damage.^{7,8}

ATM appears to be the initiating event in cells following DNA-strand breakage with the magnitude of this activation being influenced by NBS1.⁹ Both NBS1 (Nibrin) and BRCA1 (Breast cancer onset 1) migrate then to sites of DNA breaks and recruit activated ATM, which in turn phosphorylates SMC1L1.⁹ Thus, SMC1L1 phosphorylation appears to be the critical downstream event in the ATM-NBS1-BRCA1 pathway, which mediates cell survival and chromosomal stabilization after DNA damage.^{10,11} Interestingly, mutations in SMC1L1 were found in the Cornelia de Lange syndrome, a multisystem developmental disorder.¹²

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilutions should be discarded if not used within 12 hours.

Product Profile

Immunoblotting: a working dilution of 1:500-1:1,000 is recommended using HeLa nuclear extracts.

Immunoprecipitation: 5-10 μ L of the antibody immunoprecipitates SMC1L1 from HEK 293-T cell lysates.

Indirect immunofluorescence: a working dilution of 1:100-1:200 is recommended using paraformaldehyde-fixed HEK 293-T cells.

Note: In order to obtain the best results using various techniques and preparations, we recommend determining the optimal working dilutions by titration.

References

1. Losada, A., et al., *Genes & Dev.*, **19**, 1269-1287 (2006).
2. Sumara, I., et al., *J. Cell Biol.*, **151**, 749-761 (2000).
3. Hirano, T., et al., *Genes & Dev.*, **16**, 399-414 (2002).
4. Strunnikov, A.V., *J. Cell. Biol.*, **123**, 1635-1648 (1993).
5. Hagstrom, K.A., and Meyer, B.J., *Nature Rev. Genet.*, **4**, 520-534 (2003).
6. Potts, P.R., and Yu, H., *Mol. Cell. Biol.*, **25**, 7021-7032 (2005).
7. Kim, S.-T., et al., *Genes & Dev.*, **16**, 560-570 (2002).
8. Xu, Y., and Baltimore, D., *Genes & Dev.*, **10**, 2401-2410 (1996).
9. Lim, D-S., et al., *Nature*, **404**, 613-617 (2000).
10. Xu, B., et al., *Mol. Cell. Biol.*, **21**, 3445-3450 (2001).
11. Kitagawa, R., et al., *Genes Dev*, **18**, 1423-1438 (2004).
12. Mussio, A., et al., *Nat. Genet.*, **38**, 528-530 (2006).

NV,KAA,PHC 08/06-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.