


Novabiochem® coupling reagents for peptide synthesis

Novabiochem® offers a wide selection of high quality coupling reagents at competitive prices, from carbodiimide to phosphonium and aminium-based reagents. These reagents allow the efficient, fast and enantiomerization-free coupling of amino acids, even for sterically hindered residues. We reliably supply our customers with these reagents from catalog to bulk quantities.

PyOxim

Novabiochem® is very pleased to offer PyOxim.





- Cost effective alternative to COMU
- Excellent solubility in DMF and NMP
- Mediates coupling with low racemization or epimerization
- Ideal for fragment condensation and synthesis of cyclic peptides as the reagent cannot cause guanidinylation
- Likely to have low potential for causing allergic reactions
- Low or non-existent explosivity

PyOxim is used in exactly the same way as PyBOP, HBTU or HATU, but instead of generating a benzotriazolyl ester it forms the corresponding ester of Oxyma Pure

For more information or to order Novabiochem products, contact your local office

Argentina	+54 11 4546 8100	Japan	+81 0120 189 390
Australia	+61 3 9728 7600	Korea	+82 2 2185 3840
Brazil - Merck	+55 11 3346 8500	Malaysia	+6 03 7494 3688
Caribbean	See Central America	New Zealand	+64 06 356 7328
Central America	+50 2 2410-2300	Mexico	+528 18158 0600
Chile	+56 2 3400 000	Pakistan	+92 21 455 9210
China	+86 21 3222 4788	Peru	+51 1 6187 500
Colombia	+57 1 425 4770	Philippines	+63 2 815 4067
Ecuador	+593 2 2981677	Singapore	+65 6890 6660
Guatemala	See Central America	Taiwan	+886 2 2742 2788
Hong Kong	+852 2757 7569	Thailand	+66 2 667 8333
India	+91 22 6660 9184	Venezuela	+58 21 2235 1379
Indonesia	+62 21 841 3889	Vietnam	+84 8 526 0100
Israel	+972 3 938 7164		

To request a free copy of the 2010/2011 Novabiochem® Catalog or to download relevant publications, please go to www.merck4biosciences.com/literature



Product No.	Product Description	Size
851004	BOP	5 g 25 g 100 g
851054	CDI	25 g 100 g
851085	COMU	5 g 25 g 100 g
851091	DEPBT	5 g 25 g 100 g
851005	DSC	5 g 25 g 100 g
851013	HATU	5 g 25 g
851006	HBTU	5 g 25 g 100 g
851012	HCTU	5 g 25 g 100 g
851011	MSNT	1 g 5 g 25 g
851086	Oxyma Pure	25 g 100 g
851009	PyBOP	5 g 25 g 100 g
851010	PyBrOP	5 g 25 g 100 g
851087	PyCloCk	5 g 25 g 100 g
851095	PyOxim	5 g 25 g 100 g
851008	TBTU	5 g 25 g 100 g
851090	TFFH	1 g 5 g 25 g
851088	TOTU	5 g 25 g 100 g
851007	WSC	5 g 25 g