3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

Anti-acetyl-Histone H3 [Ac-Lys⁵⁶] produced in rabbit, IgG fraction of antiserum

Catalog Number SAB4200328

Product Description

Anti-acetyl-Histone H3 [Ac-Lys⁵⁶] is produced in rabbit using as immunogen a synthetic peptide containing acetylated Lys⁵⁶ of human histone H3 (GeneID: 8290), conjugated to KLH. The corresponding sequence is identical in many species including rat and mouse histone H3. Whole antiserum is purified using protein A immobilized on agarose to provide the IgG fraction of antiserum.

Anti-Acetyl-Histone H3 [Ac-Lys⁵⁶] specifically recognizes human and mouse [Ac-Lys⁵⁶] histone H3. The antibody may be used in various immunochemical techniques including immunoblotting (~15 kDa). Detection of the acetylated histone H3 [Ac-Lys⁵⁶] band by immunoblotting is specifically inhibited by the histone H3 [Ac-Lys⁵⁶] immunizing peptide, but not by the corresponding non-acetylated histone H3 peptide.

Histones are subjected to several covalent modifications, such as phosphorylation, methylation, acetylation and ubiquitination, that affect chromatin structure and regulate chromatin activity. 1,2 Histone modifications are thought to play an important role in cancer and disease.³ These modifications may alter chromatin structure and recruit downstream chromatinassociated proteins involved in transcription regulation. These in turn, may dictate dynamic transitions between transcriptionally active or silent chromatin states. Histones H3 and H4 are the predominant histones subjected to extensive covalent modifications.^{4,5} Histone H3 can be reversibly acetylated at Lys residues 9. 14. 18. 23. 27. 36 and 56. Acetylation of Lys⁵⁶ (H3-K56Ac) occurs in the N-helical region of histone H3 near the entry-exit site of the DNA superhelix, and is thought to lead to the unfolding of nucleosomes and chromatin. Acetylated Lys⁵⁶ on histone H3 has been shown to drive chromatin assembly after repair and signals for completion of repair. Histone H3 Lys56 acetylation has been linked to the core transcriptional network in human embryonic stem cells.

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide as a preservative.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilutions should be discarded if not used within 12 hours.

Product Profile

<u>Immunoblotting</u>: a working dilution of 1:1,000-1:2,000 is recommended using NIH3T3 or HeLa cell lysates.

Note: In order to obtain the best results using various techniques and preparations, we recommend determining the optimal working dilutions by titration.

References

- 1. Kouzarides, T., Cell, 128, 693-705 (2007).
- 2. Strahl, B.D., and Allis, C.D., *Nature*, **403**, 41-45 (2000).
- 3. Schneider, R., et al., *Trends Biochem.*, **27**, 396-402 (2002).
- 4. Rice, J.C., and Allis, C.D., *Curr. Opin. Cell Biol.*, **13**, 263-273 (2001).
- Garcia, B.A., et al., J. Biol. Chem., 282, 7632-7640 (2007).
- 6. Chen, C.C., et al., Cell, 134, 231-243 (2008).
- 7. Xie, W., et al., Mol. Cell, 33, 417-427 (2009).

ER,KAA,PHC 11/11-1