

QUALITY CONTROL RANGES
MILLIPLEX® MAP Human Kidney Injury Magnetic Bead Panel 2
Kit Catalog # HKI2MAG-99K
Control Catalog # HKI2-6099-2
Lot # HKI2-109 and HKI2-209



<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
EGF	Control 1	0.07 - 0.15	ng/mL
	Control 2	0.76 - 1.6	ng/mL
Lipocalin-2/NGAL	Control 1	0.17 - 0.35	ng/mL
	Control 2	1.5 - 3.1	ng/mL
Albumin	Control 1	29 – 70	ng/mL
	Control 2	266 - 553	ng/mL
Clusterin	Control 1	38 - 78	ng/mL
	Control 2	372 - 773	ng/mL
Cystatin C	Control 1	0.09 - 0.18	ng/mL
	Control 2	1.2 - 2.5	ng/mL
OPN	Control 1	5.0 - 10	ng/mL
	Control 2	57 - 119	ng/mL
α 1-Microglobulin	Control 1	6.7 - 14	ng/mL
	Control 2	71 - 148	ng/mL

**Note: The Quality Control Ranges are generated using assay buffer provided in the kit.
Quality Control Ranges in culture media are not tested.**

HKI2-6099-2: HKI2-109.209 (27 February 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA