



**QUALITY CONTROL RANGES
HUMAN ADIPONECTIN RIA KIT
CATALOG # HADP-61HK**

**STANDARD, LYOPHILIZED
CATALOG# 8061-K**

STANDARD	HUMAN ADIPONECTIN VALUE
LOT NUMBER: 18HA61	200 ng/mL

The exact value of lyophilized standard is shown above.
Reconstitute with **1 mL** distilled or deionized water.

**QUALITY CONTROLS 1 AND 2, LYOPHILIZED
CATALOG #6361-K**

Reconstitute each vial with **1 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HAD-113	6.4 – 13	ng/mL
QC 2	HAD-213	47 – 98	ng/mL

6361-K: HAD-113.213 (20 June 2024)
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA