



Specification Sheet


NovaSeptum® Accessories




NovaSeptum® accessories include a volume indicator, crimped tube caps and a holder guard. The volume indicator is used with NovaSeptum® High Purity flat bags in order to indicate the approximate volume of your sample. The crimped tube cap protects the crimped pinch pipe during transportation and storage. The holder guard protects the NovaSeptum® sterile sampling assembly throughout the process. A certificate of quality is included in each box.



Crimped Tube Cap



Holder Guard for 1.5 in.
NovAseptic® Connector



Holder Guard for 2.0 in.
NovAseptic® Connector

Specifications

Material of construction				
Volume Indicator: PVC				
Crimped Tube Caps: PVC				
Operating Temperature				
Volume Indicator:	-15 °C to 55 °C (5 °F to 131 °F)			
Crimped Tube Caps:	-20 °C to 60 °C (-4 °F to 140 °F)			
Autoclavable:	No			
Dimensions	Thickness (µm)	Volume (mL)	Length x Width (mm)	Diameter (internal)
Volume Indicator:	200 +/- 20	50	160 +/- 1 x 30 +/- 1	
	200 +/- 20	100	169 +/- 1 x 30 +/- 1	
	200 +/- 20	250	202 +/- 1 x 30 +/- 1	
	200 +/- 20	1000	325 +/- 1 x 30 +/- 1	
Crimped Tube Caps:	1.00 +/- 0.25	—	12.7 +/- 2 x 8 +/- 0.25	8.00 mm +/- 0.25 mm
Packaging				
Volume indicators are packaged in an sealed plastic bag that is packaged in a sealed outer bag and then packaged in a cardboard box.				
Ten crimped tube caps are packaged in a sealed plastic bag. Ten of these plastic bags, with a total of 100 crimped tube caps, are placed in a sealed outer bag and then packaged in a cardboard box.				

Ordering Information

Device Description	Volume (mL)	Qty/Pk	Catalogue No.
Volume Indicator	50	50	NSRULER-10050
	100	50	NSRULER-10100
	250	50	NSRULER-10250
	1000	50	NSRULER-11000
Crimped Tube Caps	--	100	NSTP-3X6

Device Description	Catalogue No.	Size	Compatible with NovaSeptum® holders	Compatible with NovAseptic® connectors
Holder Guard	NSG-NA1.5	1.5 in.	AT51/5-3D0, AT52/5-3D0	NAC-USI150, NAC-USM150, NAC-A150, NAC-I25, NAC-D32, NAC-D40
	NSG-NA2	2 in.	AT91/6-3D0, AT92/6-3D0	NAC-USI200, NAC-USM200, NAC-A200, NAC-I32, NAC-I40, NAC-D50

To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 01805 045 645

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe,
please call: +44 (0) 115 943 0840

Or visit: www.merckmillipore.com/offices

For Technical Service visit:
www.merckmillipore.com/techservice



www.merckmillipore.com