

## 67425 MDR Acinetobacter Millichrome™ plus Supplement

Selective supplement if screening is focused on MDR (multi drug resistant) *Acinetobacter*.

Package information:

5 vials (1 vial is sufficient for 1 L of final media)

Directions

- Rehydrate one vial with 5 mL of purified water.
- Add 5 mL of this solution to 1 L of the melted Acinetobacter Millichrome™ plus Agar base at 45-50 °C.
- Stir well for homogenization

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

