

## 42981 Legionella Growth Supplement (Twin Pack)

A chemical enrichment supplement recommended for enhancing growth of Legionella species.

### Composition:

(per vial, sufficient for 500 ml medium)

#### Part A:

L-Cysteine hydrochloride 200.0 mg

#### Part B:

Ferric pyrophosphate, soluble 125.0 mg

Sodium selenite 5.0 mg

Distilled water 5.0 ml

### Directions:

Rehydrate the contents of part A with 5 ml sterile distilled water. Warm up the refrigerated solution of Part B to room temperature. Aseptically add one vial of each Part A and Part B to 500 ml of sterile, rehydrated, cooled (45-50°C) Legionella Enrichment Broth, Base (Cat. No. 94423). Mix well and pour into sterile Petri plates.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

