

QUALITY CONTROL RANGES

ACTIVE GHRELIN RIA KIT CATALOG # GHRA-88HK

STANDARD, LYOPHILIZED CATALOG # 8088-K

STANDARD	ACTIVE GHRELIN VALUE
LOT NUMBER: 21A88	1959 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **2 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # 6088-K

Reconstitute each vial with **1 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	GACT-121	32 – 67	pg/mL
QC 2	GACT-221	263 – 545	pg/mL

6088-K: GACT-121.221 (rev. 12 January 2026)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.