



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sia1.com
sigma-aldrich.com

Product Information

RABBIT IgG IMMUNOPRECIPITATION REAGENT

Developed in Goat
For Immunoassay

Product No. **R 8633**
Store at 2-8 °C

Product Description

Rabbit IgG Immunoprecipitation Reagent has been specifically formulated for use in the double antibody radioimmunoassays (RIA), enzyme immunoassays (EIA), and fluoroimmunoassays (FIA) listed in the table below. The reagent is supplied ready to use for a minimum of 100 tests.

Reagents

Rabbit IgG Immunoprecipitation Reagent is a solution of goat anti-rabbit IgG antiserum in 0.01 M phosphate buffer, pH 7.4, containing 5 mM EDTA, 3.9% polyethylene glycol and 0.01% thimerosal as a preservative. Before use, gently mix by inversion.

Precautions and Disclaimer

This product is for research use only. Consult the MSDS for information regarding hazards and safe handling practices

See the individual immunoassay tracer product insert or call Sigma Technical Service for more information.

Sigma Product Number

Analyte	RIA	EIA	FIA
Acyclovir	A 7678	A 3177	N/A
Azidothymidine	A 7803	A 3302	A 9550
Chloramphenicol	C 2670	C 2545	C 1043
Chloroquine	C 1803	C 1678	N/A
Dideoxycytidine	D 0914	N/A	N/A
Deoxydihydrothymidine	D 2292	D 2167	N/A
Dideoxyinosine	D 8914	D 2417	N/A
Ubiquitin	N/A	U 9897	U 5504

JWM/VLE 4/03

Sigma brand products are sold through Sigma-Aldrich, Inc.
Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.
Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.
Please see reverse side of the invoice or packing slip.