

Cells Stained with Monoclonal Anti-a-Smooth Muscle Actin (A 2547). Staining of a-smooth muscle actin in cultured bovine smooth muscle cells using Monoclonal Anti-a Smooth Muscle Actin (A 2547, Clone No. 1A4). An intense labeling of stress fibers can be seen running in parallel in the cytoplasm and inserting at the cell membrane. The nucleus is stained red with propidium iodide. [From Dr. J. Schaper and Dr. S. Hein, Max Planck Institute, Bad Nauheim, Germany.]