

53477 CFC Selective Supplement

An antibiotic supplement for the selective isolation of *Pseudomonas* species.

Composition

(per vial, sufficient for 500 ml medium):

Cetrimide	5mg
Fucidin	5mg
Cephaloridine	25mg

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten King Agar B (Cat. No. 60786). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

