

## 17191 Park and Sanders Selective Supplement I

An antibiotic supplement for the selective isolation and enumeration of thermotolerant *Campylobacter* species.

### Composition

(per vial, sufficient for 1000 ml medium)

Vancomycin 10mg

Trimethoprim 10mg

### Directions

Rehydrate the contents of 1 vial aseptically with 5ml of 50% Ethanol. Mix well and add aseptically to 940ml of sterile, molten *Campylobacter* Selective Agar, Base (Prod No 21378) along with 50-70 ml sterile defibrinated blood. Mix well and dispense.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

