



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### MONOCLONAL ANTI-CALDESMON CLONE CALD-5 Mouse Ascites Fluid

Product No. **C 6542**

#### Product Description

Monoclonal Anti-Caldesmon (mouse IgG1 isotype) is derived from the hybridoma produced by the fusion of mouse myeloma cells and splenocytes from an immunized mouse. Turkey gizzard caldesmon was used as the immunogen. The isotype is determined using Sigma ImmunoType™ Kit (Product Code ISO-1) and by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents (Product Code ISO-2).

Monoclonal anti-Caldesmon is directed against an epitope located on the calmodulin non-binding part of the molecule. The antibody localizes the high molecular weight form of caldesmon from chicken fibroblasts when used in immunoblotting. By indirect immunofluorescence, Monoclonal Anti-Caldesmon stains frozen tissue sections from chicken, turkey, and human. Caldesmon is an actin interacting and calmodulin binding protein found in smooth muscle and other cell types. Caldesmon occurs in both a high molecular weight (120-150 kDa) and a low molecular weight (71-80 kDa) form, depending on the tissue in which it is located. Caldesmon plays a major role in the regulation of smooth muscle and non-muscle contractile events.

Monoclonal anti-Caldesmon may be used for the study of the involvement of caldesmon in contraction, cell movement, shape change, exocytosis and endocytosis.

#### Reagents

The product is provided as ascites fluid with 0.1% sodium azide as a preservative.

#### Precautions and Disclaimer

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

#### Product Profile

A working dilution of 1:100 was determined by an immunoblot assay using a chicken gizzard extract and immunoperoxidase labeling.

In order to obtain best results it is recommended that each individual user determine their optimum working dilution by titration assay.

#### Storage

For continuous use, store at 2-8 °C for up to one month. For extended storage, solution may be frozen in working aliquots. Storage in "frost-free" freezers is not recommended. Repeated freezing and thawing is not recommended. If slight turbidity occurs upon prolonged storage, clarify by centrifugation before use.

PCS/KMR 04/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.