



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-521-8956 • (314) 771-5765
Fax 800-325-5052 • (314) 771-5757
sigma-aldrich.com

Product Information

ANTI-MOUSE WHOLE SERUM

Delipidized, Whole Antiserum

Product No. **M 5774**

Product Description

Anti-Mouse Whole Serum was developed in rabbit using pooled normal mouse serum as the immunogen. The antiserum has been treated to remove lipoproteins.

Strong reactivity with normal mouse serum has been determined by immunoelectrophoresis (IEP). This antiserum has not been assayed for interspecies cross reactivity.

Reagents

Rabbit Anti-Mouse Serum is provided as a liquid containing 15 mM sodium azide as preservative.

Precautions and Disclaimer

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

Storage

For continuous use, store at 2-8 °C for up to one month. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Pcs5/00

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.