



## Technical Data Sheet

1.00246.0250

### Platin/Cobalt color reference solution (Hazen 500) Certipur®

According to DIN EN ISO 7887, ASTM D 1209-05 and DIN EN ISO 6271-1 and APHA 2120 B&C


The discoloration of water can be expressed in 'Hazen' units. The color reference solution can be used in visual comparison or calibration standards for instrumental measurements. The Hazen 500 solution is used to prepare color reference solutions with different Hazen values.

#### Benefits:

- Stable solutions with a minimum shelf life of 18 month.
- Each bottle is packed with a certificate of analysis

#### Spezifikationen:

Wavelength [nm] (Accuracy $\pm 0.2$ nm)	Absorbance (tolerance)
430	0.110 to 0.120
455	0.130 to 0.145
480	0.105 to 0.120
510	0.055 to 0.065



#### Ordering Information:

Ordering number: 1.00246.0250