## SIGMA-ALDRICH®

#### sigma-aldrich.com

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

## MARK3 (27-end), active, GST-tagged, human PRECISIO<sup>®</sup> Kinase recombinant, expressed in *Sf*9 cells

Catalog Number **M8572** Lot Number 091M0952 Storage Temperature –70 °C

Synonyms: KP78, CTAK1, PAR1A

## **Product Description**

Microtubule affinity-regulating kinase 3 (MARK3) is a member of the PAR-1/MARK kinase family, which plays critical roles in polarity and cell cycle control, and is regulated by 14-3-3 scaffolding proteins, as well as the LKB1 tumor suppressor kinase and atypical protein kinase C (PKC).<sup>1</sup> MARK3 is a dual-specificity protein kinase that controls entry into mitosis by dephosphorylating CDC2. MARK3 seems to be a positive regulator of the  $\beta$ -catenin pathway and an inhibitor of the JNK pathway. MARK3, a regulator of polarity, is also a modulator of Wnt- $\beta$ -catenin signaling, indicating a link between two important developmental pathways.<sup>2</sup>

This recombinant product was expressed by baculovirus in *Sf*9 insect cells using an N-terminal GST-tag. The gene accession number is NM 002376. It is supplied in 50 mM Tris-HCI, pH 7.5, with 150 mM NaCI, 0.25 mM DTT, 0.1 mM EGTA, 0.1 mM EDTA, 0.1 mM PMSF, and 25% glycerol.

Molecular mass: ~110 kDa

Purity: ≥70% (SDS-PAGE, see Figure 1)

Specific Activity: 614-830 nmole/min/mg (see Figure 2)

## **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

## Storage/Stability

The product ships on dry ice and storage at -70 °C is recommended. After opening, aliquot into smaller quantities and store at -70 °C. Avoid repeated handling and multiple freeze/thaw cycles.

## Figure 1.

SDS-PAGE Gel of Lot Number 091M0952: >75% (densitometry)



## Figure 2.

#### Specific Activity of Lot Number 091M0952: 722 nmole/min/mg



## Procedure

**Preparation Instructions** 

Kinase Assay Buffer – 25 mM MOPS, pH 7.2, 12.5 mM glycerol 2-phosphate, 25 mM MgCl<sub>2</sub>, 5 mM EGTA, and 2 mM EDTA. Just prior to use, add DTT to a final concentration of 0.25 mM.

Kinase Dilution Buffer – Dilute the Kinase Assay Buffer 5-fold with distilled  $H_2O$ .

Kinase Solution – Dilute the active MARK3  $(0.1 \ \mu g/\mu)$  with Kinase Dilution Buffer to the desired concentration. Note: The lot-specific specific activity plot may be used as a guideline (see Figure 2). It is recommended the researcher perform a serial dilution of active MARK3 kinase for optimal results.

10 mM ATP Stock Solution – Dissolve 55 mg of ATP in 10 ml of Kinase Assay Buffer. Store in 200  $\mu l$  aliquots at –20 °C.

 $\gamma$ -<sup>32</sup>P-ATP Assay Cocktail (250  $\mu$ M) – Combine 5.75 ml of Kinase Assay Buffer, 150  $\mu$ l of 10 mM ATP Stock Solution, 100  $\mu$ l of  $\gamma$ -<sup>32</sup>P-ATP (1 mCi/100  $\mu$ l). Store in 1 ml aliquots at –20 °C.

Substrate Solution – Dissolve the synthetic peptide substrate (KKKVSRSGLYRSPSMPENLNRPR) in water at a final concentration of 1 mg/ml.

1% phosphoric acid solution – Dilute 10 ml of concentrated phosphoric acid to a final volume of 1 L with water.

## Kinase Assay

This assay involves the use of the <sup>32</sup>P radioisotope. All institutional guidelines regarding the use of radioisotopes should be followed.

- 1. Thaw the active MARK3, Kinase Assay Buffer, Substrate Solution, and Kinase Dilution Buffer on ice. The  $\gamma$ -<sup>32</sup>P-ATP Assay Cocktail may be thawed at room temperature.
- 2. In a pre-cooled microcentrifuge tube, add the following solutions to a volume of 20 μl:
  - 10 µl of Kinase Solution
  - 5 µl of Substrate Solution
  - 5  $\mu$ l of cold water (4 °C)
- Set up a blank control as outlined in step 2, substituting 5 μl of cold water (4 °C) for the Substrate Solution.
- 4. Initiate each reaction with the addition of 5  $\mu$ l of the  $\gamma$ -<sup>32</sup>P-ATP Assay Cocktail, bringing the final reaction volume to 25  $\mu$ l. Incubate the mixture in a water bath at 30 °C for 15 minutes.
- 5. After the 15 minute incubation, stop the reaction by spotting 20  $\mu$ l of the reaction mixture onto an individually precut strip of phosphocellulose P81 paper.

- 6. Air dry the precut P81 strip and sequentially wash in the 1% phosphoric acid solution with constant gentle stirring. It is recommended the strips be washed a total of 3 times of ~10 minutes each.
- 7. Set up a radioactive control to measure the total  $\gamma^{-32}$ P-ATP counts introduced into the reaction. Spot 5 µl of the  $\gamma^{-32}$ P-ATP Assay Cocktail on a precut P81 strip. Dry the sample for 2 minutes and read the counts. Do not wash this sample.
- 8. Count the radioactivity on the P81 paper in the presence of scintillation fluid in a scintillation counter.
- 9. Determine the corrected cpm by subtracting the blank control value (see step 3) from each sample and calculate the kinase specific activity

## Calculations:

1. Specific Radioactivity (SR) of ATP (cpm/nmole)

SR =  $cpm of 5 \mu l of \gamma^{-32}P-ATP Assay Cocktail nmole of ATP$ 

cpm – value from control (step 7) nmole – 1.25 nmole (5 μl of 250 μM ATP Assay Cocktail)

2. Specific Kinase Activity (SA) (nmole/min/mg)

nmole/min/mg = 
$$\frac{\Delta \text{cpm} \times (25/20)}{\text{SR} \times \text{E} \times \text{T}}$$

SR = specific radioactivity of the ATP (cpm/nmole ATP)  $\triangle$ cpm = cpm of the sample – cpm of the blank (step 3) 25 = total reaction volume

- 20 = 101a1 reaction volume
- 20 = spot volume
- T = reaction time (minutes)
- E = amount of enzyme (mg)

## References

- Göransson, O. et al., Regulation of the polarity kinases PAR-1/MARK by 14-3-3 interaction and Phosphorylation. J. Cell Sci., **119**, 4059-4070 (2006).
- Kato, T. et al., Isolation of a novel human gene, MARKL1, homologous to MARK3 and its involvement in hepatocellular carcinogenesis. Neoplasia, 3, 4-9 (2001).

PRECISIO is a registered trademark of Sigma-Aldrich Co. LLC.

## TD,MAM 09/11-1

©2011 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.