

Product No. F-2278
Fibroblast Growth Factor-4
(FGF-4)
Human, Recombinant
Expressed in *E. coli*

Description

Fibroblast Growth Factor-4 (FGF-4), of the fibroblast growth factor family, is encoded by the K-fgf/hst oncogene. FGF-4 is also known as hst-1 (from human stomach tumor), and K-FGF (from Kaposi's sarcoma). The K-fgf/hst oncogene was originally identified by transfection into NIH3T3 cells with DNA from stomach and colon cancer,^{1,2,3} Kaposi's sarcoma^{4,5} and hepatocellular carcinoma⁶. The activation of the K-fgf oncogene results from unregulated expression of a normal gene product⁷. Unlike the two prototypes of the FGF family, FGF basic (FGF-2) and FGF acidic (FGF-1), FGF-4 is glycosylated and secreted by cells into culture medium.⁵

Fibroblast cell lines expressing FGF-4 acquire a transformed phenotype *in vitro* and are highly tumorigenic *in vivo*.⁸ FGF-4 is a potent mitogen for fibroblasts in culture and has features of a transforming growth factor.⁸ In addition, FGF-4 is mitogenic for endothelial cells and will cause morphological transformation of NIH3T3 cells.^{4,5,7} FGF-4 shares 42% sequence identity with FGF-basic (FGF-2)⁹, and both FGF-basic and FGF-4 bind to the same receptors.^{10,11} FGF-4 has a molecular weight of 22 kD.

Performance Characteristics

The biological activity of human recombinant FGF-4 is determined in a mitogenic assay by measuring ³H-thymidine incorporation in NR6-3T3 fibroblasts.¹² The EC₅₀ is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell based bioassay.

Product Information

Expressed in *E. coli*
Molecular Weight: 22 kD
Purity: ≥ 97% as determined by SDS-PAGE
EC₅₀: 0.05 - 0.3 ng/ml
Package size: 25 µg
Lyophilization Buffer: Lyophilized from a 0.2 µm-filtered solution of 20 mM sodium phosphate buffer, pH 7.0.
Carrier Protein: 1.25 mg human serum albumin (HSA).
Sterility: 0.2 µm-filtered, aseptic fill
Endotoxin: ≤ 0.1 ng/µg FGF-4

Reconstitution

Reconstitute the contents of the vial using sterile-filtered PBS containing 0.1% HSA or BSA to a concentration not less than 1 µg/ml.

Storage

Prior to reconstitution, store at -20°C. After reconstitution, store at 2-8°C for a maximum of 3 months. For extended storage, freeze in working aliquots at -70°C or -20°C. Repeated freezing and thawing is not recommended.

References

1. Sakamoto, H., et al., Proc. Natl. Acad. Sci. USA, **83**, 3997 (1986).
2. Koda, T., et al., Jpn. J. Cancer Research, **78**, 325 (1987).
3. Taria, M., et al., Proc. Natl. Acad. Sci. USA, **84**, 2980 (1987).
4. Delli-Bovi, P., et al., Proc. Natl. Acad. Sci. USA, **84**, 5660 (1987).
5. Delli-Bovi, P., et al., Cell, **50**, 729 (1987).
6. Nakagama, H., et al., Jpn. J. Cancer Research, **78**, 651 (1987).
7. Delli-Bovi, P., et al., Mol. Cell. Biol., **8**, 2933 (1988).
8. Talarico, D., et al., Molecular and Cellular Biology, **11**, 1138 (1991).
9. Bernharroch, D., et al., Isr. J. Med. Sci., **26**, 212 (1990).
10. Moscatelli, H., et al., J. Cell Biol., **109**, 2519 (1989).
11. Mansukhani, A., et al., Proc. Natl. Acad. Sci. USA, **87**, 4378 (1990).
12. Thomas, K., Methods in Enzymology, **147**, 120 (1987).

BIOHAZARD: Handle as if capable of transmitting infectious agents. Refer to MSDS.
Source material tested and found negative for antibody to HIV and HBsAG.

Sigma warrants that its products conform to the information contained in this and other Sigma publications. Purchaser must determine the suitability of the products for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale. Issued 06/96.