

## ProductInformation

### GEORG KIMMIG SELECTIVE SUPPLEMENT

Product Number **G4664**

#### Product Description

An antimicrobial supplement recommended for the selective isolation and cultivation of fungi.

#### Components

##### Item

Penicillin G	20,000 IU
Streptomycin	20.00 µg

#### Precautions and Disclaimer

For laboratory use only. Not for drug, household or other uses.

#### Preparation Instructions

Rehydrate the contents of one vial aseptically with 5 mls of sterile, distilled water. Mix well and aseptically add it to 500 mls of sterile, molten (K2635) Kimmig Fungi Agar Base. Mix well and pour into sterile petri plates.

#### Storage

Store the supplement at 2-8°C.

G4664:2/00-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.