

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Cat. No. G6633 DL-GLUTETHIMIDE--DEA SCHEDULE II

Sedative; hypnotic.

Mol. Formula: $C_{13}H_{15}NO_2$

Mol. Wt.: 217.26 (anhyd.)

m.p.: 84-86°C

CAS Registry No.: 77-21-4

Chemical Name: (±)-3-Ethyl-3-phenyl-2,6-piperidinedione

Physical Properties: White solid.

Caution: While no human toxicity data is available for this substance, it should be handled with care. Precautions

should be taken to avoid contact by all routes of exposure. RTECS No. MA4725000.

Storage: Store tightly sealed at room temperature.

Solubility: Soluble in ethanol, methanol or aqueous base; sparingly soluble in water (0.3 mg/ml). Solubility in 45% (w/v)

aqueous 2-hydroxypropyl-β-cyclodextrin (Cat. No. H-107): 6.0 mg/ml.

Disposal: Dissolve or mix the compound with a combustible solvent and burn in a chemical incinerator equipped with

an afterburner and scrubber. This product is controlled by the Drug Enforcement Administration. Appropriate security must be maintained until the substance is destroyed. Records must be kept which detail

the ultimate disposition of the material. OBSERVE ALL LOCAL, STATE AND FEDERAL LAWS.

References:

1. *Merck Index* **12th Ed.**, No. 4485.

2. Goodman and Gilman's The Pharmacological Basis of Therapeutics 9th Ed., pp. 380-383 (1996).

CONTROLLED SUBSTANCE - DEA LICENSE REQUIRED (II)