

QUALITY CONTROL RANGES
Mouse CVD Magnetic Bead Panel 2
Kit Catalog # MCVD2MAG-77K



Control Catalog # MCVDMG2-6077-2, Lot # MCVD2-106 and MCVD2-206

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>	<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
CXCL16	Control 1	175 – 365	pg/mL	sCD40L	Control 1	627 – 1302	pg/mL
	Control 2	1441 – 2992	pg/mL		Control 2	5035 – 10457	pg/mL
Endocan-1	Control 1	43 – 88	pg/mL	PIGF-2	Control 1	17 – 35	pg/mL
	Control 2	347 – 721	pg/mL		Control 2	141 – 292	pg/mL
LIGHT	Control 1	3165 – 6574	pg/mL	Troponin I	Control 1	3183 – 6611	pg/mL
	Control 2	26538 – 55118	pg/mL		Control 2	27719 – 57570	pg/mL
Follistatin	Control 1	1151 – 2391	pg/mL	Troponin T	Control 1	912 – 1894	pg/mL
	Control 2	9236 – 19182	pg/mL		Control 2	7945 – 16501	pg/mL
Oncostatin M	Control 1	153 – 318	pg/mL				
	Control 2	1227 – 2548	pg/mL				