




Product Information

SIGMA QUALITY CONTROL TEST PROCEDURE

Suitability Assay for COLLAGEN

PRINCIPLE:



CONDITIONS: T = 37°C, pH = 7.4, A_{570nm}, Light path = 1 cm

METHOD: Colorimetric

REAGENTS:

- A. 50 mM TES Buffer with 0.36 mM Calcium Chloride, pH 7.4 at 37°C
(Prepare 1000 ml in deionized water using TES Free Acid, Sigma Prod. No. T-1375, and Calcium Chloride, Dihydrate, Sigma Prod. No. C-3881. Adjust the pH to 7.4 at 37°C with 1 M NaOH.)
- B. Collagen
(Different types of Collagen will produce varying amounts of enzyme activity when used as a substrate for collagenase.)
- C. Collagenase Enzyme Solution
(Immediately before use, prepare a solution containing 0.1 mg/ml Collagenase, Sigma Prod. No. C-0130, in Buffer A.)
- D. Ethylene Glycol Monoethyl Ether
(Use Ethylene Glycol Monoethyl Ether, Sigma Prod. No. E-2632.)
- E. 4% (w/v) Ninhydrin Solution
(Prepare 100 ml in Reagent D, using Ninhydrin, Sigma Prod. No. N-4876.)
- F. 200 mM Citrate Buffer with 0.16% (w/v) Stannous Chloride, pH 5.0 at 25°C
(Prepare 100 ml in deionized water using Citric Acid, Free Acid, Anhydrous, Sigma Prod. No. C-0759. Adjust to pH 5.0 at 25°C with 1 M NaOH. Then add the Stannous Chloride, Anhydrous, Sigma Prod. No. S-2752.)

Suitability Assay for COLLAGEN

REAGENTS: (continued)

- G. 50% (v/v) 1-Propanol Solution
(Prepare 100 ml in deionized water using 1-Propanol, Sigma Stock No. 29,328-8.)
- H. Ninhydrin Color Reagent (NCR)
(Immediately before use, combine equal volumes of Reagent E and Reagent F.)
- I. 10 mM Hydrochloric Acid Solution
(Prepare 50 ml in deionized water using Hydrochloric Acid, Sigma Prod. H-7020.)
- K. 4.0 mM L-Leucine Standard Solution (Std Soln)
(Prepare 20 ml in Reagent I using L-Leucine, Sigma Prod. No. L-8000. **PREPARE FRESH.**)

PROCEDURE:

Weigh the following reagent into suitable containers:

	Test	Blank
Reagent B (Collagen)	25.00 mg	25.00 mg

Then add (in milliliters) the following reagent:

Reagent A (Buffer)	5.00	5.00
--------------------	------	------

Incubate the vials at 37°C until equilibrated. Then add:

Reagent A (Buffer)	-----	0.10
Reagent C (Enzyme Solution)	0.10	-----

Mix well and incubate at 37°C. Swirl the containers for 10 - 15 seconds at 1.5 and 3.5 hours. After 5 hours, filter the contents of the containers through a Whatman #54 filter paper or a 0.8 µm syringe filter into clean containers. Use the filtrates for color development.

Suitability Assay for COLLAGEN

COLOR DEVELOPMENT:

Standard Curve:

Prepare a standard curve by pipetting the following reagents (in milliliters) into suitable vials.

		<u>Std 1</u>	<u>Std 2</u>	<u>Std 3</u>	<u>Std 4</u>	<u>Std Blank</u>
Reagent K (Std Soln)	0.05	0.10	0.15	0.20	0.00	
Deionized Water		0.15	0.10	0.05	0.00	0.20
Reagent H (NCR)		2.00	2.00	2.00	2.00	2.00

Sample:

Pipette (in milliliters) the following reagents into suitable vials:

	<u>Test</u>	<u>Blank</u>
Test Filtrate	0.20	-----
Blank Filtrate	-----	0.20
Reagent H (NCR)	2.00	2.00

Mix well and place vented caps on each vial. Place the vials in a boiling water bath for 20 minutes. Remove the vials and allow to cool to room temperature. Add 10 ml of Reagent G (50% 1-Propanol) to each vial. Mix well and transfer the vial contents to suitable cuvettes. Determine the absorbance at 570 nm for each of the vials using a suitable spectrophotometer.

CALCULATIONS:

Standard Curve:

$$\Delta A_{570\text{nm}} \text{ Standard} = A_{570\text{nm}} \text{ Standard} - A_{570\text{nm}} \text{ Standard blank}$$

Prepare a standard curve by plotting the $\Delta A_{570\text{nm}}$ of the L-Leucine Standard Solution versus micromoles of L-Leucine.

Sample Determination:

$$\Delta A_{570\text{nm}} \text{ Sample} = A_{570\text{nm}} \text{ Test} - A_{570\text{nm}} \text{ Sample blank}$$

Suitability Assay for COLLAGEN

CALCULATIONS: (continued)

Determine the μ moles of L-Leucine equivalents liberated using the Standard curve.

$$\text{Units/ml enzyme} = \frac{(\mu\text{moles of L-Leucine equivalents liberated}) (5.1) (df)}{(0.2) (0.1)}$$

df = Dilution factor

5.1 = Total volume (in milliliters) of Assay

0.2 = Total volume (in milliliter) of sample used in
Colorimetric Determination

0.1 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

UNIT DEFINITION:

One unit liberates peptides from collagen equivalent in ninhydrin color to 1.0 μ mole of leucine in 5 hours at pH 7.4 at 37°C in the presence of calcium ions.

SPECIFICATION:

Suitable for use as a substrate for Collagenase.

FINAL ASSAY CONCENTRATION:

In a 5.10 ml reaction mix, the final concentrations are
50 mM TES, 0.36 mM calcium chloride, 25 mg collagen and
0.01 mg collagenase.

REFERENCES:

Moore, S. and Stein, W.H. (1948) *J. Biol. Chem.* **176**, 367-388

Mandl, I., MacLennan, J.D., Howes, E.L., DeBellis, R.H., and Sohler, A. (1953) *Journal of Clinical Investigation* **32**,
1323-1329

Suitability Assay for COLLAGEN

NOTES:

1. This assay is based on the cited references.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.