



## QUALITY CONTROL RANGES

### HUMAN HIGH MOLECULAR WEIGHT (HMW) ADIPONECTIN ELISA KIT CATALOG #: EZHMWAN-65K

#### STANDARD, LYOPHILIZED CATALOG #: E8065-K

STANDARD	High Molecular Weight Adiponectin Value
LOT NUMBER: E865-4K	200 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.25mL** distilled or deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #: E6065-K

Reconstitute each vial with **0.25mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HMW-104	6.9 – 14	ng/mL
QC 2	HMW-204	29 – 61	ng/mL

E6065-K: HMW-104.204 (10 November 2022)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)  
The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA