

Product No. D-4908

Lot 096H4850

Anti-Dog Whole Serum

Developed in Rabbit
Delipidized, Whole Antiserum

Anti-Dog Whole Serum was developed in rabbit using pooled normal dog serum as the immunogen. The antiserum has been treated to remove lipoproteins. Rabbit Anti-Dog Serum is provided as a liquid containing 0.1% sodium azide (see MSDS)* as preservative.

Specificity

Strong reactivity with normal dog serum has been determined by immunoelectrophoresis (IEP). This antiserum has not been assayed for interspecies cross reactivity.

Protein Concentration = 61.8 mg/ml by Biuret.

Storage

For continuous use, store at 2-8°C. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

*Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.