

Product Information

Hoechst Stain Solution

1.2 g/L BisBenzimide H 33258

CATALOG NO. 59412C

Description

Hoechst Stain Solution is prepared with 1.2 g/L BisBenzimide H33258 fluorochrome stain in cell culture grade water and sterile filtered for use in mycoplasma staining procedures.

Precautions

Use aseptic technique when handling this stain. This product is for research or for further manufacturing use. THIS PRODUCT IS NOT INTENDED FOR HUMAN OR THERAPEUTIC USE.

Storage

Solution should be stored at 2 to 8 C protected from light.

Indications of Deterioration

Solution should be clear and free of particulate and flocculent material. Do not use if solution is cloudy or contains precipitate. Other evidence of deterioration may include color change, pH shift and degradation of physical or performance characteristics. Do not use the product after the expiration date.

Preparation Instructions

Supplements, such as antibiotics, can be added to the sterilized solution using aseptic technique. Storage conditions and shelf life of the supplemented product may be affected by the nature of the supplements.

Characteristics

Appearance

Clear solution

pH (as supplied)

Refer to Certificate of Analysis

Sterility

No microbial growth detected

United States

SAFC Biosciences, Inc.
13804 W. 107th Street
Lenexa, Kansas 66215
USA
Phone +1 913-469-5580
Toll free-USA 1 800-255-6032
Fax +1 913-469-5584
E-mail info-na@sial.com

Europe

SAFC Biosciences Ltd.
Smeaton Road, West Portway
Andover, Hampshire SP10 3LF
UNITED KINGDOM
Phone +44 (0)1264-333311
Fax +44 (0)1264-332412
E-mail info-eu@sial.com

Asia Pacific

SAFC Biosciences Pty. Ltd.
18-20 Export Drive
Brooklyn, Victoria 3025
AUSTRALIA
Phone +61 (0)3-9362-4500
Toll free-AUS 1 800-200-404
Fax +61 (0)3-9315-1656
E-mail info-ap@sial.com

Warranty, Limitation of Remedies

S AFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC Biosciences expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

© 2006 SAFC Biosciences, Inc.

Issued September 2006 P59412
0805 1005 0406

United States

SAFC Biosciences, Inc.
13804 W. 107th Street
Lenexa, Kansas 66215
USA

Phone	+1 913-469-5580
Toll free-USA	1 800-255-6032
Fax	+1 913-469-5584
E-mail	info-na@sial.com

Europe

SAFC Biosciences Ltd.
Smeaton Road, West Portway
Andover, Hampshire SP10 3LF
UNITED KINGDOM

Phone	+44 (0)1264-333311
Fax	+44 (0)1264-332412
E-mail	info-eu@sial.com

Asia Pacific

SAFC Biosciences Pty. Ltd.
18-20 Export Drive
Brooklyn, Victoria 3025
AUSTRALIA

Phone	+61 (0)3-9362-4500
Toll free-AUS	1 800-200-404
Fax	+61 (0)3-9315-1656
E-mail	info-ap@sial.com

www.safcbiosciences.com