




- **EZBio[®]pure TR Media Bottles** are available in volumes; 125mL, 250mL, 500mL, and 1L
- **NEW** ergonomic bottle shape with molded finger grips
- GL45 VersaCap[®] made of polypropylene
- **EZBio[®]pure TR Media Bottles** available in either PS or PETG materials
- “Pure” refers to type of resin used to manufacture the bottle; USP Class VI resins
- Manufactured in a Class 100,000 clean-room
- Manufactured in an ISO 13485 Facility in North America
- **Gamma Sterilized in USA** 



Foxx Life Sciences’ **EZBio[®]pure Titanium Round (TR) Media Bottles**, have an enhanced the media bottle design. To stabilize the media bottles we lowered the bottle height, angled shoulder for smooth and complete pouring of media, and molded finger grips to improve the ability to grasp the bottle even with gloved hands. Easy to read raised graduation marks help in dispensing media.


The **EZBio[®]pure TR Media Bottles** are made of either, Polystyrene (PS) or Polyethylene terephthalate glycol (PETG). Both materials use USP Class VI resins and are crystal clear.

Each bottle has a GL45 neck that is compatible with our patented GL45 VersaCap[®] assemblies as well as other standard GL45 caps. The GL45 VersaCap[®] allows for added media security and enables easy handling of the bottle even with gloved hands. Additionally, the innovative cap features a flat writing surface, easy for labeling.

EZBio[®]pure TR Media Bottles are also designed to handle full vacuum and can be used with **EZBio[®]pure Feed Systems**.

- 1 GL45 VersaCap[®] Technology
- 2 Forward-facing, large, raised graduation marks for easy reading
- 3 Molded finger grips to improve the ability to grasp the bottle even with gloved hands





	Item Number	Volume	VersaCap [®] Assembly	Sterile or Non-Sterile	U/M	Description
Bottle Only (PETG)	507-1631-FLS	125mL	No Cap	Non-Sterile	48/CS	EZBio [®] pure TR Media Bottle, PETG, 125mL, GL45 Neck w/ No Cap
	507-1831-FLS	250mL	No Cap	Non-Sterile	48/CS	EZBio [®] pure TR Media Bottle, PETG, 250mL, GL45 Neck w/ No Cap
	507-1931-FLS	500mL	No Cap	Non-Sterile	24/CS	EZBio [®] pure TR Media Bottle, PETG, 500mL, GL45 Neck w/ No Cap
	507-1A31-FLS	1L	No Cap	Non-Sterile	24/CS	EZBio [®] pure TR Media Bottle, PETG, 1L, GL45 Neck w/ No Cap
	517-1632-FLS	125mL	Cap Only	Sterile	48/CS	EZBio [®] pure TR Media Bottle, PETG, 125mL, GL45 Neck w/ Cap
	517-1832-FLS	250mL	Cap Only	Sterile	48/CS	EZBio [®] pure TR Media Bottle, PETG, 250mL, GL45 Neck w/ Cap
	517-1932-FLS	500mL	Cap Only	Sterile	24/CS	EZBio [®] pure TR Media Bottle, PETG, 500mL, GL45 Neck w/ Cap
	517-1A32-FLS	1L	Cap Only	Sterile	24/CS	EZBio [®] pure TR Media Bottle, PETG, 1L, GL45 Neck w/ Cap
Bottle Only (PS)	506-1631-FLS	125mL	No Cap	Non-Sterile	48/CS	EZBio [®] pure TR Media Bottle, PS, 125mL, GL45 Neck w/ No Cap
	506-1831-FLS	250mL	No Cap	Non-Sterile	48/CS	EZBio [®] pure TR Media Bottle, PS, 250mL, GL45 Neck w/ No Cap
	506-1931-FLS	500mL	No Cap	Non-Sterile	24/CS	EZBio [®] pure TR Media Bottle, PS, 500mL, GL45 Neck w/ No Cap
	506-1A31-FLS	1L	No Cap	Non-Sterile	24/CS	EZBio [®] pure TR Media Bottle, PS, 1000mL, GL45 Neck w/ No Cap
	516-1632-FLS	125mL	Cap Only	Sterile	48/CS	EZBio [®] pure TR Media Bottle, PS, 125mL, GL45 Neck w/ Cap
	516-1832-FLS	250mL	Cap Only	Sterile	48/CS	EZBio [®] pure TR Media Bottle, PS, 250mL, GL45 Neck w/ Cap
	516-1932-FLS	500mL	Cap Only	Sterile	24/CS	EZBio [®] pure TR Media Bottle, PS, 500mL, GL45 Neck w/ Cap
	516-1A32-FLS	1L	Cap Only	Sterile	24/CS	EZBio [®] pure TR Media Bottle, PS, 1000mL, GL45 Neck w/ Cap
Bottle Assemblies	527-4400-FLS	125mL	Cap Assembly	Sterile	10/CS	125mL Bottle Assy., PETG, Plug and Vented w/ Diptube
	527-5400-FLS	250mL	Cap Assembly	Sterile	10/CS	250mL Bottle Assy., PETG, Plug and Vented w/ Diptube
	527-6400-FLS	500mL	Cap Assembly	Sterile	10/CS	500mL Bottle Assy., PETG, Plug and Vented w/ Diptube
	527-7400-FLS	1L	Cap Assembly	Sterile	10/CS	1L Bottle Assy., PETG, Plug and Vented w/ Diptube
	526-4400-FLS	125mL	Cap Assembly	Sterile	10/CS	125mL Bottle Assy., PS, Plug and Vented w/ Diptube
	526-5400-FLS	250mL	Cap Assembly	Sterile	10/CS	250mL Bottle Assy., PS, Plug and Vented w/ Diptube
	526-6400-FLS	500mL	Cap Assembly	Sterile	10/CS	500mL Bottle Assy., PS, Plug and Vented w/ Diptube
	526-7400-FLS	1L	Cap Assembly	Sterile	10/CS	1L Bottle Assy., PS, Plug and Vented w/ Diptube
Cap Assemblies	520-4400-FLS	125mL	Cap Assembly	Sterile	20/CS	125mL GL45 Cap Assy., PP, Plug and Vented w/ Diptube
	520-5400-FLS	250mL	Cap Assembly	Sterile	20/CS	250mL GL45 Cap Assy., PP, Plug and Vented w/ Diptube
	520-6400-FLS	500mL	Cap Assembly	Sterile	20/CS	500mL GL45 Cap Assy., PP, Plug and Vented w/ Diptube
	520-7400-FLS	1L	Cap Assembly	Sterile	20/CS	1L Cap GL45 Assy., PP, Plug and Vented w/ Diptube