



## RABBIT ANTI-RAT CYTOCHROME P450 CYP4A1, CYP4A2 and CYP4A3 POLYCLONAL ANTIBODY

<b>CATALOG NUMBER:</b>	AB1280
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	100 µL
<b>SPECIFICITY:</b>	Binds specifically to cytochrome P450 CYP4A1, CYP4A2 and CYP4A3 in rat hepatic microsomal fraction.
<b>IMMUNOGEN:</b>	Synthetic peptide
<b>APPLICATIONS:</b>	Western blot Immunohistochemistry Optimal working dilutions must be determined by end user.
<b>FORMAT:</b>	Rabbit antiserum
<b>PRESENTATION:</b>	Liquid, no preservatives.
<b>STORAGE/HANDLING:</b>	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

*For research use only; not for use as a diagnostic*