

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

## **ProductInformation**

## ASPIRIN SIGMA REFERENCE STANDARD

Product No. A3160

Storage Temperature: room

Synonyms: Benzoic acid, 2-(acetyloxy)-.

Salicylic acid acetate.

Formula: C<sub>9</sub>H<sub>8</sub>O<sub>4</sub>

Formula Weight: 180.16

## **Product Description**

This product meets USP specifications (Current through USP 27, 2nd Supplement), and is traceable to USP Reference Standard lot H.


Dry over silica gel for 5 hours before using. Keep container tightly closed.

TEST	USP 27 SPECIFICATIONS	Result Sigma R.S. Lot 084K8802
Identification A	a violet-red color is produced	passes
Identification B Infrared absorption spectrum KBr dispersion (197K)	compares to standard	conforms to USP RS lot H
Loss on drying (731)	not more than 0.5%	passes
Readily carbonizable substances (271)	no more color than <i>Matching Fluid Q</i>	passes
Residue on ignition (281)	not more than 0.05%	passes
Substances insoluble in sodium carbonate TS	clear solution	Passes
Chloride (221)	not more than 0.014%	Passes
Sulfate	not more than 0.04%	Passes
Heavy Metals (231)	not more than 10 μg per g	Passes
Limit of free salicylic acid	not more than 0.1%	Passes
Organic volatile impurities  Method IV (467)	meets requirements	Meets requirements
Assay	not less than 99.5% and not more than 100.5% calculated on the dried basis	100.0%

See reverse side for IR spectrum of Sigma Reference Standard.

## **ASPIRIN** SIGMA REFERENCE STANDARD

Product No. A3160 Lot 84K8802



INFRARED SPECTRUM

Wavenumbers (cm<sup>-1</sup>)